


No.	Plan	Article	Qty
1	E00066473	2998022	1
2	Y003413	1463254	1
3	Y007779	1607954	1
4	Y003414	1463354	1
5	E00052508	2998018	1
6	Lever M8x15	1006661	2
7	Tube cap 15x15	1005899	1
8	M6x20	1006999	4

* Use:
WEICONLOCK AN 302-43

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

PLEGO L55 HEADREST				2998012
				Plan No. E00066564/A1.1
Scale: 1:4	Drafted	Verified	Confirmed	
Date	18.08.2023	18.08.2023	18.08.2023	
Initials	BK	WM	TK	